

# BVI

## 40 Years of Innovation and Quality Cannulas

Supporting All Areas of Ophthalmic Surgery

- A market leader of ophthalmic cannulas for over 40 years
- Full range of single-use cannulas developed to support all ophthalmic surgical procedures
- A history of innovation to meet your diverse surgical needs



Visitec®  
Cannulas

# Irrigation

## Anterior Chamber Cannula

Angled 45 degrees, 6 mm from end.




581284

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	581642
● 27G x 7/8 in, 0.4 x 22 mm	581284
● 30G x 7/8 in, 0.3 x 22 mm	581285

## Viscoflow Cannula

Angled 45 Degrees, 9 mm from end.



585009

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585009
● 25G x 7/8 in, 0.5 x 22 mm	585045
● 27G x 7/8 in, 0.4 x 22 mm	585026
● 30G x 7/8 in, 0.3 x 22 mm	585046

## Anterior Chamber (Rycroft)


Angled 45 degrees, 4 mm from end.



585025

Qty: 10/Bx  
Brand: Visitec®


Description	SKU #
● 21G x 7/8 in, 0.8 x 22 mm	581281
● 23G x 7/8 in, 0.6 x 22 mm	585025
● 25G x 7/8 in, 0.5 x 22 mm	581639
● 27G x 7/8 in, 0.4 x 22 mm	581280
● 30G x 7/8 in, 0.3 x 22 mm	581273



# Refractive

## LASIK Cannula

Intended to irrigate and/or manipulate corneal tissue. 25 Degrees, formed 7 mm from end.



585278

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585278

## LASIK Irrigating Cannula (Solomon)

Intended to irrigate and/or manipulate corneal tissue. Thin wall cannula, angled tip features three 0.25 mm ports on 8 mm extension.




585280

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585280

## LASIK Irrigating Cannula (Lindstrom)

Intended to irrigate and/or manipulate corneal tissue. Features three 0.25 mm ports located on 7 mm extension with an 11 mm radius, 45 degrees.




585250

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585250

## LASIK Flap Spatula/Irrigator (Hersh)


Intended to irrigate and/or manipulate corneal tissue. 45 Degrees, 9 mm from spatulated tip.



585252

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	585252



# Anesthesia

## Retrobulbar (Atkinson)

Intended for injecting anaesthetic agents into the muscle cone.




581638

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 1-1/2 in, 0.6 x 38 mm	581638
● 25G x 1-1/2 in, 0.5 x 38 mm	581637

## Retrobulbar/Peribulbar (Atkinson)



585108

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 1-1/4 in, 0.6 x 32 mm	585001
● 25G x 1-1/4 in, 0.5 x 32 mm	585108
● 27G x 1-1/4 in, 0.4 x 32 mm	585109

## Retrobulbar (Sharp Point)

Intended for injecting anaesthetic agents into the muscle cone.




581275

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 1-1/2 in, 0.5 x 38 mm	581275

## Sub-Tenon's Anesthesia Cannula


Intended for administering anaesthetic agents into the posterior sub-tenon's space.



585176

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 19G x 1 in, 1.1 x 25 mm	585176



ISO Hub Gauge Chart

19G	20G	21G	22G	23G	24G	25G	26G	27G	30G
-----	-----	-----	-----	-----	-----	-----	-----	-----	-----

# Hydrodissection

## Cortical Cleaving Hydrodissector (Flat End)

45 degrees, 11 mm from end. For cleaving cortex from capsular bag.



585158

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585158

## Micro J Shape Hydrodissector (Pearce)

Hook is 2.25 mm long x 1.75 mm wide.




585029

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585029

## Hydrodissection Cannula (Akahoshi Style)

Tapered tip, Thin wall, Angled 1 mm from end.



585198

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585198

## Nucleus Hydrodissector

Tip is flattened on the horizontal plane.




585155

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G 7/8 in, 0.4 x 22 mm	585155

## Micro J Shape Cannula

Hook is 1.5 mm long x 1.75 mm wide. 0.25 mm open end port.




585003

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585003

## Nucleus Hydrodissector (Buratto)

Formed, 1 mm bent tip.



585199

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585199



### Nucleus Hydrodissector (Chang)

90 Degree bend, 1 mm flattened tip.



585294

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 33 degrees	585294
● 27G x 7/10 in, 0.4 x 18 mm, Bevel angle 47 degrees	585295

### Nucleus Hydrodissector (Pearce)

Angled 35 degrees, 8 mm from end.




585036

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585037
● 27G x 7/8 in, 0.4 x 22 mm with 8 mm bend	585036

### Nucleus Hydrodissector (Sauter)

Tip is flattened on the vertical plane.



585099

Qty: 10/Bx  
Brand: Visitec®


Description	SKU #
● 27G x 7/8 in, 0.4 x 22 mm	585099



# Capsule Polishers

## Capsule Polisher (Kratz)

0.3 mm Top port, 35 degrees, 9 mm from end.




581605

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 7/8 in, 0.6 x 22 mm	581605

## Flexible Tip Cannula

Angled 45 degrees, 7 mm from end, 10 mm sleeve.




585017

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 7/8 in, 0.5 x 22 mm	585017
● 27G x 7/8 in, 0.4 x 22 mm	585005

## Olive Tip Polisher

0.33 mm end port, 30 degrees, 10 mm from end and micro-etched 2 mm from end. Olive shape tip is 2 mm long x 1 mm wide.



585143

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 23G x 1 in, 0.6 x 25 mm	585143

# Cystotomes

## Irrigating Cystotome

Formed cystotomes are bent to conform with the convexity of the crystalline lens.




581610

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 5/8 in, 0.5 x 16 mm	581610
● 27G x 5/8 in, 0.4 x 16 mm	581618

## Irrigating Cystotome

Formed 8.5 mm from bend to tip. Shortened curve designed for deep set eyes or small pupils.



585008

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G 5/8 in, 0.5 x 16 mm	585004
● 27G 5/8 in, 0.4 x 16 mm	585008

## Irrigating Cystotome (Blumenthal)

Formed 8.5 mm from bend to tip.



585042

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 27G x 5/8 in, 0.4 x 16 mm	585042

# Cortex Extractors

## Cortex Extractor

0.3 mm end port angled 10 degrees. Flat end opening has rounded edges.



585106

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 24G x 1/2 in, 0.55 x 12.5 mm	585106

## Cortex Extractor

0.4 mm top port.



585193

Qty: 10/Bx  
Brand: Visitec®


Description	SKU #
● 23G x 5/8 in, 0.6 x 16 mm	585193



# Lacrimal

## Lacrimal Cannula

Intended to irrigate and probe the lacrimal system. Curved.




581615

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 21G x 1-3/8 in, 0.8 x 35 mm	581615
● 26G x 1-1/8 in, 0.45 x 28 mm	581276

## Lacrimal Cannula

Straight, 13 mm Length.



585068

Qty: 10/Bx  
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in, 0.5 x 13 mm	585068



FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Questions: [bvimedical.com/customer-support](http://bvimedical.com/customer-support)



BVI and all other trademarks (unless noted otherwise) are property of BVI. BVI ©2022

1536846-02 IN