









Mini-Blade and Handle Portfolio

Innovative Solutions for Specialty Surgery









For over 75 years, BVI has been the industry leader for meeting the needs of specialty surgeons for a broad range of surgical procedures. Each blade is color coded for easy identification.

Grindless Beaver Mini-Blades*		
Description		SKU #
Grindless, Rounded tip, Sharp on one side, Green, 20/Bx		Beaver6400
Grindless, Lamellar, Angled 60 degrees, Smooth undersurface, Sharp all around, Bevel up, Black, 3/Bx		Beaver6600
Grindless, Lamellar, Straight, Smooth undersurface, Sharp all around, Black, 20/Bx		Beaver6610
Grindless, Smaller alternative to #15 general blade, Orange, 20/Bx		Beaver6700
Grindless, Straight, Sharp all around, 180 Degree cutting surface, Double bevel blade, Blue, 20/Bx		Beaver6900
Grindless, Cleft palate blade, Angled 60 degrees, Sharp all around, 180 Degree cutting surface, Double bevel blade, Blue, 6/Bx		Beaver6910


* Patent D636082. Additional US and International patents pending.

Beaver® Mini-Blades

Beaver® Mini-Blades can be used with the following Beaver handles, no collet required

Beaver® Handles		
Description		SKU #
Round, Knurled, Stainless, 13 cm L x .5 cm D		371350
Round, Knurled, Stainless, 13 cm L x .66 cm D		371360
Round, 21 cm L x .4 cm D, No chuck collet included		371395
Round, Knurled, Stainless, 10 cm L x .5 cm D		374310
Round, Knurled, Stainless, 5 cm L x .5 cm D		374311
Round, Knurled, Stainless, 15.5 cm L x .5 cm D		374320
Hex, Stainless, 10 cm L x .5 cm D		374330
Hex, Stainless, 15.5 cm L x .5 cm D		374340

Mini-Blades Packaging



FOR COMPLETE PRODUCT INSTRUCTIONS, PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Tel: 866-906-8080 Fax: 866-906-4304
 Email: customersupport@bvimedical.com



BVI, BVI Logo and all other trademarks are property of BVI © 2019

1536847-01