

BVI

Providing Relief for Patients with Dry Eye or Lacrimal Disorders

Parasol[®] | Micro Flow[™]



**Permanent
Parasol[®]**



**Partial
Micro Flow[™]**

Single-Use Punctal Dilators and Probes

Visitec® Punctal Dilator and Probe, Khan

Dual function instrument allows for single handed transition from dilating to probing without release of lateral lid traction.



QTY: 3/Bx
Brand: Visitec®

Description	SKU #
25 mm/1 in medium taper, Blunt, Wilder style dilator mated with a size 0 Bowman lacrimal probe, 0.7 mm x 45 mm/22G x 1 3/4 in	585124

Punctum Dilator, Fine/Medium

0.4 mm x 21.5 mm Tip,
0.6 mm x 18 mm Tip,
77 mm Handle




Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Punctum Dilator, Fine/Medium	MMSU1181S

Punctum Dilator/Seeker

0.45 mm x 21.5 mm Tip,
5 mm Pointed tip,
18 mm Straight shaft, 77 mm Handle



Qty: 20/Bx
Brand: Malosa®

Description	SKU #
Punctum Dilator/Seeker	MMSU1193S



Permanent Occlusion

Permanent occlusion is indicated for chronic dry eye conditions. All occluder products are made of medical grade silicone and are reversible.

Insertion

The hollow nose collapses on insertion, eliminating the need for punctal dilation in most cases.


Straight Forward Sizing

Two sizes, small and medium, fit virtually all patients. X-Small and large sizes are available for customization.

Retention

Parasol is designed for excellent retention, as demonstrated by an independent clinical study.¹

Parasol® Size Comparison Chart



Parasol® Pre-Loaded Punctal Occluders Sterile

Occluders are packaged in a convenient 1 pair tray, are sterile and pre-loaded on a single-use inserter.



A14-201




Parasol® Pre-Loaded Punctal Occluders Sterile	SKU #
X-Small, 0.25 - 0.35 mm	A14-200
Small, 0.35 - 0.65 mm	A14-201
Medium, 0.6 - 0.85 mm	A14-202
Large, 0.9 mm+	A14-203

¹ C. McCabe, MD, PhD, FACS Plugs Reduce Dry-Eye Symptoms, Improve Vision. November 2009, Review of Ophthalmology.

Parasol® Quickload Punctal Occluders Sterile

Occluders are sterile and pre-mounted on a stand for Quickload convenience.

QTY: 16/Bx



Parasol® Quickload Punctal Occluders Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201QL
Medium, 0.6 - 0.85 mm	A14-202QL
Large, 0.9 mm+	A14-203QL
Combo small and medium	A14-205QL

Parasol® Bulk Punctal Occluders Non-Sterile

Occluders are non-sterile and packaged 2 occluders per poly bag.

QTY: 20/Bx



Parasol® Bulk Punctal Occluders Non-Sterile	SKU #
Small, 0.35 - 0.65 mm	A14-201BN
Medium, 0.6 - 0.85 mm	A14-202BN
Combo small and medium	A14-209BN

Partial Occlusion

Micro Flow™ Punctal Occluders

Designed to provide a straight-forward, reversible means of partial occlusion. Micro Flow™ Punctal Occluders support lacrimal drainage care. Suggest to use when total occlusion is not ideal. Occluders are packaged in a convenient 1 pair tray, are sterile and pre-loaded on a single-use inserter.



A14-401



QTY: 2/Bx

Micro Flow™ Pre-Loaded Punctal Occluders Sterile	SKU #
Sterile, Small, 0.4 - 0.55 mm	A14-401
Sterile, Medium, 0.55 - 0.7 mm	A14-402
Sterile, Large, 0.7 - 0.85 mm	A14-403

QTY: 10/Bx

Micro Flow™ Punctal Occluders Non-sterile	SKU #
Non-Sterile, Bulk, Small, 0.4 - 0.55 mm	A20-401BN
Non-Sterile, Bulk, Medium, 0.55 - 0.7 mm	A20-402BN
Non-Sterile, Bulk, Large, 0.7 - 0.85 mm	A20-403BN



Other Related Occlusion Options

Curved Tying Forceps

0.25 mm Tips, 5 mm Curved tying platforms,
80 mm Handle



Qty: 20/Bx
Brand: Malosa®

Description	SKU #	SKU # Clips
Curved Tying Forceps	MMSU1210S	MMSU1210CS


Lacrimal Cannula



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 25G x 1/2 in/0.5 x 13 mm	585068

Lacrimal Cannula



Qty: 10/Bx
Brand: Visitec®

Description	SKU #
● 26G x 1 1/8 in/0.45 x 28 mm	581276



FOR COMPLETE PRODUCT INSTRUCTIONS,
PLEASE REFER TO THE BVI IFU OR PRODUCT INSERT

Questions: bvimedical.com/customer-support



BVI
bvimedical.com